INSTRUCTIONS FOR USE

Stride Duo - OA KNEE BRACE

art no 70656

Provider information

Range of motion adjustments:

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
EN The Stride Duo comes with standard 5° extension stops. Additional flexion and extension stops may be installed if they are prescribed by your physician or needed for optimal performance or fit of the brace. These should be inserted by the technician who fits your brace. Optional 0°, 10°, 15°, 20°, and 25° extension stops and optional 45°, 60°, 75°, and 90° flexion stops are included with each brace.

IMPORTANT: Both hinges must be set to same extension and/or flexion settings. For 0° extension limiting, 0° extension stops must be used.

TO CHANGE AN EXTENSION STOP: Unscrew the screw on the front of the hinge cap. Shake out the 5° extension stop. Insert another extension stop into the front of the hinge and replace the screw. Repeat for the opposite hinge.

TO INSERT A FLEXION STOP: Insert a flexion stop into the back of the hinge and use the extra screw that is attached to the flexion/extension stop tree to secure. Repeat for the opposite hinge.

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EN Successively increase the time of using the brace, start using the brace for 30 minutes the first day. To obtain best possible effect, wear the brace during activities that strain the knee, such as walking and standing. Arthritis caused pain may take a few weeks to improve.

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EN For single patient use only. If a stinging sensation, rash or other changes in skin color occur, discontinue use and contact your doctor or other medically trained person. Serious incidents shall be reported to NordiCare and the competent authority.

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EN Mild to severe unicompartmental knee osteoarthritis with ligament instabilities. Increased rehabilitative and functional support for diagnoses that require reduced medial or lateral joint loads. **Contraindications:** Not intended for people with weight >110 kg.

For more information please consult our product catalogue or visit www.nordicare.se. For further questions, contact your doctor, health practitioner or medical equivalent - or NordiCare: +46 42-35 22 20, info@nordicare.se.

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EN Salt water, sand, dirt and other elements can cause damage. Thoroughly rinse brace and air dry. Liners should be removed and hand washed as needed using cold water and mild detergent, rinse thoroughly and air dry. Lubricate hinges with dry silicone spray as needed. Periodically check the condition of the Velcro to ensure they are free from grease and damage.

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EN Aluminum bar, neoprene knee bag. Lined with suede-like material of polyester. The silicone on the inside for anti-slip.

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EN Stride Duo has double uprights for increased knee stability. Used for pain relief and to facilitate and improve knee function in case of osteoarthritis.

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EN 1. Unfasten all four posterior (behind leg) straps and fold strap ends back onto themselves. In a seated position with leg slightly bent, place brace on leg so that center of hinge is lined up with top of knee cap (1 & 2). Hinge center is marked by the gray extension stop.

2. Attach straps in sequence indicated on strap ends. Hold brace at hinge to prevent movement while tightening straps.

• **Strap 1:** Fasten upper calf strap and position above calf muscle; just below the crease behind the knee (3).

• **Strap 2:** Fasten lower calf strap.

• **Strap 3:** Fasten lower thigh strap.

• **Strap 4:** Fasten upper thigh strap securely but not too tightly. Overtightening may restrict circulation and/or result in downward migration of brace.

• **Strap 5:** Fasten front strap (5) and tighten securely (4).

NOTE: From side view, hinges should sit slightly posterior (behind) to midline of knee. Adjustment of straps 1, 3, and 5 will influence this alignment.

3. When positioning the Stride Duo off-the-shelf brace, it should first be adjusted to neutral. This is achieved by setting the brace angle so the cuffs and hinge rest on leg without gapping or pressure. Using the correction key, adjust Q-hinges on both sides of the brace until they are aligned with the patient's leg. If the patient's leg has a varus or valgus deformity, the brace should be angled at the same degree. There should be no pressure on either side of the knee. For patients with severe deformity of the knee, this step may need to be completed before tightening all straps. Stand up and take a few steps. Brace will settle slightly to align with center of kneecap. If adjustment is needed, hold brace to prevent rotation while straps are repositioned. If the brace drops below knee center, reapply the brace.

OA CORRECTION: Initial application and adjustment should be done by your healthcare practitioner or orthotic professional. Re-adjustment by the patient after initial fitting is rarely necessary.

Apply correction for medial OA:

• Use the correction key and turn the Q-hinge on the inside of the knee 1/4 turn (90°) **CLOCKWISE** to open the orthosis away from the leg (6a).

• Turn the Q-hinge on the outside of the knee 1/4 turn **COUNTERCLOCKWISE** to apply correction to the knee (6b).

Apply correction for lateral OA:

• Use the correction key and turn the Q-hinge on the outside of the knee 1/4 turn (90°) **CLOCKWISE** to open the orthosis away from the leg (7a).

• Turn the Q-hinge on the inside of the knee 1/4 turn **COUNTERCLOCKWISE** to apply correction to the knee (7b).


Additional adjustments of medial/lateral OA:

• Let the patient walk around a little. They should feel a slight pressure on the outside/inside of the knee. If they do not feel a light pressure, correct a little more by turning the Q-hinge on the outside/inside another 1/8 turn (45°). If they feel an uncomfortable pressure, reduce the Q-hinge on the inside by 1/8 turn. Do not overcorrect. Do not make corrections when the patient stands up.


ATC ADJUSTMENTS:

NOTE: Applying correction to the Q-hinges may cause buckling or straining to one or more straps of the active thigh cuff. If there is a wrinkle in the ATC, tighten that strap one notch. If there is more tension in one strap relative to the other, loosen that strap one notch. Trim the end of the ATC so that it does not stick out outside the brace frame. Stride Duo braces with the ATC may be adjusted further if the uprights are not positioned along the midline of the thigh or if there is gapping between the soft tissue and the ATC straps. ATC adjustment should be done by your healthcare professional or brace technician. Remove the screws for each strap. Adjust as needed to fit flush against the thigh with the side flanges sitting mid-medial and mid-lateral or slightly posterior. Replace screws. Trim as needed.


ATC MODELS: If desired, install ATC pivot clip where the ATC components cross (8). The pivot clip allows the ATC to pivot in the center during use and is there as an option if the crossing straps are moving too much or shift out of their proper place on the leg. Stride Duo comes with two types of ATC straps – one is more flexible, one is more rigid. Use the type recommended for you.

APPLY CORRECTION FOR MEDIAL OA

APPLY CORRECTION FOR LATERAL OA

ATC PIVOT CLIP