INFORMATION FOR

Cirrus - AIRWALKER

art no 80215

Indications

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EN If an increased height is wanted for the plantarflexion a heel wedge can be used as an accessory. Five independent wedges are assembled together to the desired height. For ökad komfort finns strumpor i 2-pack som tillbehör. Hygiene protection sold in 10-pack.

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EN Immobilization of foot and ankle. For the treatment of soft tissue injuries, sprains, stable fractures and an achilles tendon rupture. Even for postoperative stabilization of the ankle.

Contraindications: Not intended for diabetic patients. Not intended for unstable fracture. **Risks:** No known.
Size XS and S, maximum weight 90 kg
Size M, L and XL maximum weight 135 kg

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Material Composition

62% Polypropylene, 5% Rubber, 11% Nylon, 8% Polyurethane (PU) foam, 3% Metal, 6% Thermoplastic Polyurethane (TPU), 3% Ethylene vinyl acetate (EVA), 2% Paper Pulp. Latex-free.

For more information, see productleaflet or www.nordicare.se. You can also contact your fitter or Nordicare: +46 42-35 22 20 / DK 7070-1907, info@nordicare.se.

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EN Always read all instructions carefully before using this product. When in doubt, consult your doctor/device fitter. It is important that you follow the general and individual user instructions you receive from your doctor/device fitter throughout your rehabilitation period.

Prescription, initial fitting and recommendation on how to use this product should be made by medically trained person. If changes in color of skin or other blood circulation related problems and skin irritation occurs, immediately remove the brace. The brace is only for one patient use. Serious incidents shall be reported to NordiCare and the competent authority.

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EN The shell is cleaned with a damp cloth. The padding should be hand washed and dried flat. Do not tumble dry.

MD Class I **CE**

This packaging can be recycled as paper and the bag as plastic.

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






EN Instruction Fitting

- When applying Cirrus for first time open all straps, **A** remove anterior shell **B** then remove liner **C**. Toe guard **D** can be adjusted or removed, see step **6** below.
- Apply liner ensuring the heel is pushed fully into the liner. Fasten hook closures on foot and leg sections ensuring they are snug.
- Place foot/liner into walker ensuring heel is against posterior border of walker. Apply pressure to fix liner to shell.
- Using hook attachments fix front plastic shell to front of liner.
- Fasten straps, from toes up, ensuring they are firm but not too tight.
- Toe guard can be adjusted or removed if required.
- For narrow legs or if increased compression is need, inserting both ends of ankle strap through slots in walker section as shown.
- Remove the transport guards over the pumps. Press the pumps to inflate air cells to the required level. Ensure patient feedback is used to determine comfort. Both pumps can be pressed simultaneously for balanced pressure.
- Open the cap and press the valve to deflate.

Integrated pumps on both sides regulate compressions independently of each other.

The circular bands can be slipped through the shell to allow for maximum compression.

Removable and adjustable toe protection.

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After Initial Fitting

Remove Cirrus:

- Open all straps and remove anterior plastic shell.
- Open liner and remove lower leg.
- Liner should remain attached to walker.

Apply Cirrus:

- Insert lower leg into the walker assembly ensuring heel is well back.
- Fasten liner ensuring it is snug.
- Using hook attachments fix shell to front of liner.
- Fasten straps, from toes up, ensuring they are firm but not too tight.