

Alex – facial bandage

NEW



NordiCare[®] +

NORDICARE ORTOPEDI + REHAB AB • Solrosvägen 1 • SE-263 62 Viken, Sverige
Tel: 042-35 22 20 • Fax: 042-35 22 21 • info@nordicare.se • www.nordicare.se

Alex – facial bandage

NEW

A simple and flexible bandage for fixation and compression after facial surgery

- Effective compression after for example facelift, liposuction or ear surgery
- Adjustable fit and compression with velcro that attaches over the entire bandage. This means the bandage fits many different sizes and shapes.
- Two supplied narrower straps can be used for extra fixation.
- The material breathes and is soft and comfortable for the skin.
- The fine elastic helps the bandage to follow the contours of the head.



The bandage can be placed around the chin/neck with the narrower straps at the back or front for extra fixation.



The bandage can also be used without the narrow straps.



For compression of ears or hairline, Alex can be placed over the forehead and around the back of the head, with or without narrow bands.

Alex facial bandage

Size	Width	Colour	Art.no
Universal	7,5 cm	white	51108-09-50

CE

Latex free

